



PCR (Polymerase Chain Reaction)

www.dauids-bio.com (Custom Antibodies)

www.dauids-science.de (Lab Material)

- 1 - Introduction

PCR (Polymerase Chain Reaction) is a fundamental technique employed to amplify a specific segment of DNA. It was developed in 1983 and has since become one of the most utilized methodologies in the fields of molecular biology and genetics. By allowing the exponential amplification of specific DNA fragments *in vitro*, PCR enables researchers to generate abundant quantities of genetic material for further analysis. This enzymatic assay relies on the use of two flanking oligonucleotides called primers, which selectively bind to targeted regions of the DNA template, as well as a thermostable DNA polymerase.

The process of PCR involves cycling through different temperatures. Initially, the double-stranded template DNA is denatured and separated into its single-stranded strains. During the annealing step, the primers bind to the complementary single-stranded DNA molecules. In the extension step, the primers are elongated and DNA synthesis is performed by the DNA polymerase.

PCR finds extensive applications in research, particularly in the study of gene expression, mutations and genetic variations. It is an invaluable tool for detecting the presence of DNA from pathogens such as viruses and bacteria.

- 2 - Material

| Material | |
|------------|---|
| Template | DNA of interest with the region for amplification |
| Master Mix | Distilled nucleases free water (nuclease free) |
| | DNA Polymerase (thermostable) |
| | Forward primer |
| | Reverse primer |
| | dNTP Mix (dATP, dCTP, dGTP, dTTP) |
| | Polymerase Buffer |
| Vial | Thin walled 200 µl tubes |

- 3 - Method

Sample Preparation

- Gather the necessary number of PCR tubes
- Prepare the PCR master mix with enough volume for the desired number of reactions
Please refer to "Sample Approach" for required components
- Distribute the master mix evenly among sterile PCR tubes
- Add the DNA template to each tube
- Thoroughly mix the PCR reaction solution by gentle pipetting
- Briefly centrifuge the tubes
- Place the tubes into a thermal cycler and proceed with the amplification program

Sample Approach

| Reagents | For 25 µl Approach | Final Concentration |
|-------------------------|--------------------|---------------------|
| 10x Polymerase Buffer | 2.5 µl | 1x |
| 10 mM dNTP Mix | 0.5 µl | 200 µM |
| Thermostable Polymerase | 0.25 µl | 2.5 Units |
| 10 µM Forward Primer | 0.5 µl | 0.2 – 0.4 µM |
| 10 µM Reverse Primer | 0.5 µl | 0.2 – 0.4 µM |
| Sample DNA | X µl | 5 ng |
| Water | ad. 25 µl | |

PCR Program

- Initiate the standard amplification program on the thermal cycler
- Set the "Primer Annealing" temperature to the melting temperature minus 5°C [TM -5] of utilized the primer
Common primer annealing temperatures are between 48 – 60 °C
- Set the "Fragment Elongation" duration to expected fragment length
Common duration for unknown fragments is 1 minute

Thermo Cycler Program

| Steps | Temperature | Duration |
|---------------------------|-------------|-----------|
| Initial Denaturation | 95°C | 3 min |
| <i>Start Cycle of 35x</i> | | |
| Template Denaturation | 95°C | 30 sec |
| Primer Annealing | [TM – 5] °C | 30 sec |
| Fragment Elongation | 72°C | 30 sec/kb |
| <i>End Cycle</i> | | |

| | | |
|------------------|------|-------|
| Final Elongation | 72°C | 5 min |
| Storage | 4°C | - |

Detection by Agarose gel

- Weigh the appropriate amount of agarose
Concentration of 1.5% is suitable for PCR products ranging from 100 bp to 2000 bp
- Dissolve the agarose in TAE buffer, pH 8
- Heat the agarose solution to completely dissolve the agarose
Please refer to the safety datasheets and handle Ethidium bromide with care
- Allow the agarose solution to cool down slightly
- Add 0.5 µg/ml Ethidium bromide
- *Approximately 1 µl per 100 ml of agarose solution*
- Carefully load your DNA samples into the gel
- Limit the current to a maximum of 5 mA for a single gel
- Run the gel at 100 V for 40 minutes
- To visualize the DNA bands, use a UV lamp at a wavelength of 302-365 nm
Take necessary precautions, as exposure to UV light can be harmful to the skin and eyes

- 4 - Trouble Shooting

| Reason | Solution |
|----------------------------|--|
| No Amplification | <p>Ensure good DNA template quality and concentration.</p> <p>Verify primer design and sequences are appropriate.</p> <p>Confirm the thermal cycler settings (temperatures, ramp rates, and cycling program).</p> <p>Optimize annealing temperature by performing a temperature gradient or adjusting primer concentrations.</p> |
| Non-specific Bands | <p>Optimize primer design to ensure specificity</p> <p>Consider a Hot-Start PCR to minimize non-specific amplification, by increasing the initial denaturation to 15 minutes</p> |
| High Background | <p>Check for contamination in reagents or equipment</p> <p>Utilize high-quality DNA template</p> <p>Use nuclease-free water</p> <p>Reduce the cycle number</p> |
| GC rich templates | <p>Optimize annealing temperature</p> <p>Consider using additives such as DMSO</p> <p>Design new primers with adjusted melting temperatures</p> <p>Use a specialized DNA Polymerase with enhanced GC-rich amplification capabilities</p> |
| Primer Dimer Contamination | <p>Optimize primer design to reduce complementarity</p> <p>Adjust primer concentrations and optimize annealing temperature</p> |