

# WesternBlot

www.davids-bio.com

## SOLUTIONS

| SDS Sample Buffer   | Running Buffer  | Transfer Buffer (wet)   | Blocking Buffer  |
|---|---|---|--|
| 125 mM Tris-HCl<br>4% w/v SDS<br>20% Glycerol<br>10% β-Mercaptoethanol<br>0.004 % Bromphenol Blue<br>pH 6.8 | 25 mM Tris base<br>200 mM Glycine<br>0.1% SDS<br>pH 8.5 | 25 mM Tris base<br>200 mM Glycine<br>15% Methanol<br>0.1% SDS<br>pH 8.5 | 1 x TBS<br>0.1% Tween-20<br>5% BSA <b>or</b> Milk Powder<br><br>Stirr and filtrate |
| Wash Buffer   |   |   |  |
| 1 x TBS<br>0.1% Tween-20  |   |   |  |

## SAMPLE PREPARATION

- Y Use an appropriate amount of sample solution by determining the protein concentration by Bradford or 280 nm.
- Y Dilute the sample with “SDS Sample Buffer” 1:2 to 1:5
- Y Heat the sample for 95°C for 5 min
- Y Store the samples at -20°C or continue with the protocol

## SDS-GEL

- Y 10% gels can be used for most proteins in the range of 10 – 100 kDa
- Y Pipet equal amounts of protein (10 – 200 ng) or cell lysates (10 – 40 µg) to the gel
- Y You may want to use pre-stained markers to have a control for membrane blotting
- Y Run the gel: 80 – 140 V for 1 – 2 hours. You may want to start with 80 V and turn it up after 10 – 20 min

## ELECTROTRANSFER TO A MEMBRANE

- Y Use Nitrocellulose or PVDF. PVDF must be wetted with Methanol for 1 min before usage. Remove the Methanol with Transfer Buffer.
- Y You may want to dye the membrane with Ponceau red to check if the transfer is successful
- Y Follow the instructions of your Blotter or Membrane for optimal results



## STAINING WITH ANTIBODIES

- Y **Blocking:** 1 h at room temperature with blocking solution (Alternatively block at 2 – 8°C over night).
- Y **Wash:** Membrane 3 times with Wash Buffer
- Y **Primary Antibody:** Incubate the membrane with your primary antibody over night at 2 – 8°C  
*Affinity purified antibodies: 0.1 – 20 µg/ml (i. e. 1:50 from a 0.5 mg/ml antibody solution)*  
*ProteinA purified antibodies: 1 – 200 µg/ml (i. e. 1:100 from a 10 mg/ml antibody solution)*
- Y **Wash:** Membrane 3 times with Wash Buffer
- Y **Secondary Antibody:** Incubate 1 h at room temperature with the recommended dilution of the antibody (refer to manual of the secondary antibody)
- Y **Wash:** Membrane 3 times with Wash Buffer
- Y **Substrate:** Follow the instructions of your secondary antibody or substrate (i. e. TMB, pNPP)